# **Q800R3** DNA Shearing Sonicator

# QSONICA

- High throughput
- Uses standard sample tubes
- 2 year warranty

## **APPLICATIONS:**

- DNA shearing for Next-gen sequencing
- Chromatin shearing
- ChIP
- ChIP-seq
- RNA-seq
- Protein extraction
- Cell lysis

### **KEY INSTALLATIONS:**

- Broad Institute
- NIH
- Harvard Medical School
- Memorial Sloan Kettering
- Max Planck Institute
- Ludwig Institute UCSD
- University of Cambridge
- Washington Univ. St. Louis
- California Institute of Technology
- Shanghai Inst. for Biological Sciences

Qsonica has over 45 years of experience manufacturing ultrasonic equipment with thousands of customers around the world. All products are designed and manufactured in the USA and include a full 2-year warranty.

GEONICA

GEONICA

The Q800R3 is our 2nd generation DNA and Chromatin shearing system. Improvements include an enhanced user-friendly design and quieter operation while remaining thousands of dollars less than the competition.

Samples can be sheared to a range of fragment sizes (150bp – 3kb). The system is compatible with commercially available sample tubes and is capable of processing up to 18 samples at one time. Multiple tube rack options can accommodate 50ul - 1ml sample volumes using a variety of standard sample tubes.

The Q800R3 is a complete package including a chiller to control temperature, digital operating system with programmable memory, high intensity ultrasonics and sound reducing enclosure.

# Chromatin Shearing

Example protocols and results are based on customer feedback.

### Sample Protocols:

### Mammalian Chromatin

Cell Type: HEK 293T (2) 10cm dishes 70-80% confluent

Total Sample Volume: 300ul

Fixation Time: 1% Formaldehyde, 13 min

Sonicator Amplitude Setting: 70%

Sonication Pulse Rate: 15 seconds On, 45 seconds Off

Total Sonication On Time: 30 min

#### Yeast Chromatin

Cell Type: Wild type S. pombe cells grown to an OD of 1.3-1.5 in YEA.

Total Sample Volume: 300ul

Fixation Time: 1% Formaldehyde, 15 min

Sonicator Amplitude Setting: 100%

Sonication Pulse Rate: 20 seconds On, 40 seconds Off

Total Sonication On Time: 30 min

### C. elegans Chromatin

Cell Type: C. elegans culture; Nuclear extract 1mg/ml

Total Sample Volume: 500ul

Fixation Time: 2% Formaldehyde, 20 min

Sonicator Amplitude Setting: 70%

Sonication Pulse Rate: 30 seconds On, 30 seconds Off

Total Sonication On Time: 20 min

### Sample Results:



### Mammalian Chromatin Prep & Chip



Chr10 Primer sequences:Forward - TCCTTCTCCCAACAATCAGCReverse - GATGTCGCTCCGAATCTTG Anitbodies Used: H3K9me3 (abcamab8898), H3K9Ac (upstate 07-352), H3K36me3(abcam ab9050) H3 (abccamab1791)

\* Average of two independent replicate chromatins in pane



## Yeast Chromatin Prep & Chip



\* Note: Lane M is the NEB 1kb Plus ladder

Sample protocols and publications for additional species and cell lines are available in the literature section of www.sonicator.com

# **DNA Shearing**

Example protocols and results are based on customer feedback.

# Sample Protocols:

### **Bacterial Genomic DNA**

Cell Type / Concentration: E.coli / 250ng Sonication Pulse Rate: 15 seconds On, 15 seconds Off Total Sample Volume: 200ul in 0.5ml tubes Total Sonication On Time: As indicated below Sonicator Amplitude Setting: 20%

#### FU ±20 100 00 60 40 12 m 6 min 4 min 20 min 20 100 150 200 35 300 400 500 600 1000 2000 10360 (bo)

## Ful 100 00 60 40 12 min 20 (bp) 100 150 400 500 600 ÷. 300 1000 2000 10390

# Human Skeletal Myoblast Genomic DNA

Cell Type / Concentration: LHCN-M2 Skeletal Myoblast / 200ng Sonication Pulse Rate: 15 seconds On, 15 seconds Off Total Sample Volume: 200ul in 0.5ml tubes Total Sonication On Time: As indicated below **Sonicator Amplitude Setting:** 20%







# Fragmentation Over Time

### Testing Independent Replicates

www.sonicator.com/dna



### Sample Tube Rack Options for the Q800R3 System

#4256 8 tube holder (1.5mL Polystyrene tubes)



#4255 12 tube holder (0.5mL PCR tubes)





### Technical Specifications:

Power Rating:	750 Watts			
Frequency:	20 kHz			
Programmable Timer:	1 second to 10 hours			
Voltage:	110V*, 50/60Hz			
*Specify desired voltage for export				

#### **Dimensions:**

Generator	8.00" W	Х	15.25" L	Х	8.50" H
Enclosure	11.50" W	Х	12.00" L	Х	20.00" H
Chiller	11.00" W	Х	13.00" L	Х	13.00" H

Additional information on our website:

- Sample Protocols
- Optimization tips
- Video
- Customer Feedback
- Recent Publications



Qsonica, LLC. 53 Church Hill Road, Newtown, CT 06470 Phone: 203.426.0101 E-mail: Info@Sonicator.com

www.sonicator.com/dna

